

Advancing the Arizona State University Knowledge Enterprise

Inventors

David Brafman Xiao Wang Nicholas Brookhouser Stefan Tekel Kylie Standage-Beier

Contact

Jovan Heusser jovan.heusser@skysonginnovat ions.com

Base-Edited Isogenic hPSC Line Generation (BIG-TREE)

Human pluripotent stem cells (hPSCs) are a vital tool not only in research and disease modeling but also in developmental biology and regenerative medicine because of their ability to self-renew and differentiate. Current genome editing tactics in hPSCs require the incorporation of deleterious double-stranded DNA breaks followed by inefficient homology-directed repair (HDR). These methods can result in DNA insertion or deletion issues, apoptosis, oncogenic mutations and other negative outcomes.

Researchers at Arizona State University have developed a novel method using CRISPR tools to co-target genomic loci and an episomal reporter to enable singlenucleotide genomic changes in hPSCs without HDR. This method, called baseedited isogenic hPSC line generation using a transient reporter for editing enrichment (BIG-TREE), enables the rapid generation of clonal isogenic hPSC lines. It produces single-nucleotide editing efficiencies greater than 80% across multiple hPSC lines and allows for the efficient generation of loss-of-function hPSC lines and efficient multiplex editing of hPSCs at multiple independent loci.

BIG-TREE advances the implementation of base-editing technologies in hPSCs for use in developmental biology, disease modeling, drug screening, cell-based therapies and more.

Potential Applications

- Isogenic cell line engineering
- Generating knockout cell lines
- Developmental biology
- Disease modeling
- Drug screening
- Regenerative medicine

- Greater than 80% efficiencies across multiple hPSC lines
- Allows for the precise and efficient base editing of hPSCs
- Efficient generation of loss-of-function hPSC lines
- Efficient editing of multiple loci and across several independent hPSC lines
- Bulk enrichment of base-edited cell populations including hPSCs

• Fast and efficient generation of clonal isogenic hPSC lines with homozygous and heterozygous single base pair edits

- Rapid engineering of isogenic hPSC lines
- o Establishes in vitro models to assess pathogenic risk/disease mechanisms

• Allows for biallelic/multiplexed targeting without the need for sequential retargeting

Compatible with off the shelf chemical transfection reagents

o Doesn't require cloning of complex viral constructs/use of specialized cell transfection systems

• Can be used in conjunction with other base editing variants that have altered PAM specificities and editing windows

For more information about this opportunity, please see

Brookhouseer et al - Stem Cell Reports - 2020

For more information about the inventor(s) and their research, please see

Dr. Brafman's departmental webpage

Dr. Wang's laboratory webpage